



COMPONENTS



PRODUCT CODES



AM0440

Antegrade IR Platform without Extension x 1



AM0410

IR Platform Extension x 1




AM0430

Antegrade IR Platform + Extension x1

UNBOXING

1



AM0440


Antegrade IR Platform Only

AM0410

Extension Only

AM0430

Antegrade IR Platform + Extension




CLEANING



Disinfect and wipe clean before/after each use.
DO NOT USE any other sterilisation methods.

ASSEMBLY



For an online training visit:
AdeptMedicalTraining.com



For service and repairs visit
AdeptMedical.com/repairs



For Approved Disinfectants PDF visit:
AdeptMedical.com/antegrade-ir-platform

MANUFACTURER

Adept Medical Ltd
2-6 McDonald St, Morningside
Auckland, New Zealand
www.adeptmedical.co.nz




EC REP

EC REPRESENTATIVE

Emergo Europe
Prinsessegracht 20
2514 AP The Hague
The Netherlands




SETUP



DRAPESUPPORT (AM1000)* | ARMSURE (AM0600)*

Optional



*Drape Support and ArmSure are not included with the Antegrade IR Platform.

